

Anti-Gldc Rabbit Polymab® For Wb, Ihc-P, If/Icc, Elisa - P23378

Item Number: CM0009874



Introduction

High-performance Anti-GLDC Rabbit PolymAb® validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Targets glycine dehydrogenase, mitochondrial, for metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GLDC Rabbit PolymAb®
Remarks/Alias	GCE; GCSP; HYG1
Species	Human
Gene ID	2731
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-749 of human GLDC (NP_000161.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23378
Protein Weight	113kDa
Shipping	Ice bag