

Anti-Glucosylceramidase Beta (Gba) Monoclonal Antibody For Wb, Ihc-P, Elisa - P04062

Item Number: CM0002024



Introduction

Rabbit monoclonal antibody targeting human Glucosylceramidase beta (GBA) validated for Western blot, IHC paraffin, and ELISA. Cross-reacts with human and rat. Ideal for research on lysosomal lipid metabolism and Gaucher disease.

[Learn More](#)

Parameter	Specification
Product Name	Glucosylceramidase beta (GBA) Rabbit mAb
Remarks/Alias	GBA; GCB; GLUC; Glucosylceramidase beta (GBA)
Species	Human
Gene ID (Human)	2629
Gene ID	2629
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 437-536 of human Glucosylceramidase beta (GBA) (P04062)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P04062
Protein Weight	60kDa
Shipping	Ice bag