


Anti-GLUD1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P00367

Item Number: CM0028247



Introduction

Rabbit polyclonal antibody targeting human GLUD1, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Suitable for studies on glutamate dehydrogenase, mitochondrial metabolism, and insulin homeostasis in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	GLUD1 Rabbit pAb
Remarks/Alias	GDH; GDH1; GLUD; hGDH1; GLUD1
Species	Human
Gene ID (Human)	2746
Gene ID	2746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 54-242 of human GLUD1 (NP_005262.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00367
Protein Weight	61kDa
Shipping	Ice bag