

Anti-Glut2/Slc2A2 Rabbit Monoclonal Antibody For Wb And Elisa - P11168

Item Number: CM0003672

Introduction

This rabbit monoclonal antibody targets GLUT2/SLC2A2, a facilitative glucose transporter. Validated for WB and ELISA, it cross-reacts with human, mouse, and rat samples, making it ideal for diabetes and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GLUT2/SLC2A2 Rabbit mAb
Remarks/Alias	GLUT2; GLUT2/SLC2A2
Species	Human
GeneID (Human)	6514
GeneID	6514
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GLUT2/SLC2A2 (P11168).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11168
Protein Weight	57kDa
Shipping	Ice bag