


Anti-Glycophorin C (Gypc) Monoclonal Antibody For Wb, Ihc-P, Elisa - P04921

Item Number: CM0002289



Introduction

Rabbit monoclonal antibody targeting Glycophorin C (GYPC) for Western blot, immunohistochemistry, and ELISA. Shows cross-reactivity with human, mouse, and rat. Useful for red cell membrane studies and blood group antigen research.

[Learn More](#)

Parameter	Value
Product Name	Glycophorin C (GYPC) Rabbit mAb
Remarks/Alias	GE; GPC; GPD; GYPD; CD236; PAS-2; CD236R; PAS-2'; Glycophorin C (GYPC)
Species	Human
GeneID (Human)	2995
GeneID	2995
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 49-128 of human Glycophorin C (GYPC) (P04921).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04921
Protein Weight	14kDa
Shipping	Ice bag