


Anti-Gnai2 Rabbit Polyclonal Antibody For Wb, Elisa - P04899

Item Number: CM0025335



Introduction

Anti-GNAI2 rabbit polyclonal antibody for Western blot and ELISA, reactive against human and mouse GNAI2. Mediates GPCR signaling, inhibits cAMP production, and plays roles in angiogenesis and inflammation.

[Learn More](#)

Parameter	Value
Product Name	GNAI2 Rabbit pAb
Remarks/Alias	GIP; HG1C; GNAI2B; H_LUCA15.1; H_LUCA16.1; GNAI2
Species	Human
Gene ID (Human)	2771
Gene ID	2771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human GNAI2 (NP_002061.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04899
Protein Weight	40kDa
Shipping	Ice bag