

Anti-Gnaq Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P50148

Item Number: CM0002553

Introduction

Rabbit monoclonal antibody targeting human GNAQ (G-alpha-q), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for GPCR signaling research and diagnostic raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	GNAQ Rabbit mAb
Remarks/Alias	GAQ; SWS; CMAL; CMC1; G-ALPHA-q; GNAQ
Species	Human
GeneID (Human)	2776
GeneID	2776
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human GNAQ (NP_002063.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50148
Protein Weight	42kDa
Shipping	Ice bag