

Anti-Gng10 Rabbit Polyclonal Antibody For Wb, Elisa - P50151

Item Number: CM0027407

Introduction

High-quality anti-GNG10 rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human, mouse, and rat GNG10, a 7 kDa G protein subunit involved in transmembrane signaling.

[Learn More](#)

Parameter	Value
Product Name	GNG10 Rabbit pAb
Remarks/Alias	GNG10
Species	Human
GenelD (human)	2790
GenelD	2790
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human GNG10 (NP_001017998.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50151
Protein Weight	7kDa
Shipping	Ice bag