


Anti-Gosr2 Rabbit Polyclonal Antibody For Wb, Elisa - O14653

Item Number: CM0026298



Introduction

High-quality rabbit polyclonal antibody against human GOSR2 (membrin, GS27), validated for WB and ELISA. Suitable for Golgi transport research and SNARE complex studies.

[Learn More](#)

Parameter	Value
Product Name	GOSR2 Rabbit pAb
Remarks/Alias	Bos1; EPM6; GS27; MYOS; GOSR2
Species	Human
GeneID (Human)	9570
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14653
Protein Weight	25kDa
Shipping	Ice bag