


Anti-Got1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P17174

Item Number: CM0002542



Introduction

High-performance rabbit monoclonal antibody against human GOT1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat; 46 kDa protein. Ideal for neuroscience and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	GOT1 Rabbit mAb
Remarks/Alias	AST; AST1; SGOT; cCAT; GIG18; cAspAT; ASTQTL1; GOT1
Species	Human
Gene ID	2805
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human GOT1 (NP_002070.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17174
Protein Weight	46kDa
Shipping	Ice bag