

Anti-Gp6/Gpvi Monoclonal Antibody For Fc, Elisa - Q9Hcn6

Item Number: CM0004079



Introduction

Rabbit monoclonal antibody targeting human GP6/GPVI, validated for flow cytometry and ELISA. Suitable for studying platelet collagen receptor function in thrombosis and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	GP6/GPVI Rabbit mAb
Remarks/Alias	GPVI; GPVI; BDPLT11; GP6/GPVI
Species	Human
GeneID (Human)	51206
GeneID	51206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-267 of human GP6/GPVI (NP_057447.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	Q9HCN6
Protein Weight	37kDa
Shipping	Ice bag