

Anti-Gpc4 Rabbit Polyclonal Antibody For Wb And Elisa - 075487

Item Number: CM0016575

Introduction

This rabbit polyclonal antibody targets GPC4 (K-glypican), suitable for WB and ELISA applications. Cross-reacts with human and mouse species. Protein weight ~62kDa. Useful for investigating glypican-4 function in tubule development and central nervous system.

[Learn More](#)

Parameter	Description
Product Name	GPC4 Rabbit pAb
Remarks/Alias	KPTS; K-glypican; GPC4
Species	Human
GeneID (human)	2239
GeneID	2239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 357-556 of human GPC4 (NP_001439.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075487
Protein Weight	62kDa
Shipping	Ice bag