

Anti-Gper1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99527

Item Number: CM0017714



Introduction

Rabbit polyclonal GPER1 antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat GPER1. Suitable for studying estrogen signaling, cardiovascular and reproductive research. Sw-Prot Q99527.

[Learn More](#)

Parameter	Value
Product Name	GPER1 Rabbit pAb
Remarks/Alias	mER; CEPR; GPER; DRY12; FEG-1; GPR30; LERGU; LyGPR; CMKRL2; LERGU2; GPCR-Br
Species	Human
GeneID (Human)	2852
GeneID	2852
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 307-375 of human GPER1 (NP_001496.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99527
Protein Weight	42kDa
Shipping	Ice bag