

Anti-GpIX/CD42A/Gp9 Rabbit Monoclonal Antibody For Wb, IHC-P, Elisa - P14770

Item Number: CM0006582



Introduction

This rabbit monoclonal antibody targets human GPIX/CD42a (GP9), a subunit of the GPIb-V-IX complex, and is validated for WB, IHC-P, ELISA. Ideal for hemostasis and platelet adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	GPIX/CD42a/GP9 Rabbit mAb
Remarks/Alias	GPIX; CD42a; GPIX/CD42a/GP9
Species	Human
GeneID	2815
Immunogen	Recombinant protein of human GPIX/CD42a/GP9
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P14770
Protein Weight	19kDa
Shipping	Ice bag