

Anti-Gpnmb Rabbit Polyclonal Antibody For Wb And Elisa - Q14956

Item Number: CM0018545

Introduction

Rabbit polyclonal antibody against human GPNMB (NMB, HGFIN), validated for western blot and ELISA. Shows cross-reactivity with mouse and rat. Detects endogenous protein at ~64 kDa, useful for oncology and dermatology studies.

[Learn More](#)

Parameter	Value
Product Name	GPNMB Rabbit pAb
Remarks/Alias	NMB; HGFIN; PLCA3; GPNMB
Species	Human
Gene ID (Human)	10457
Gene ID	10457
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-486 of human GPNMB (NP_001005340.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14956
Protein Weight	64kDa
Shipping	Ice bag