

Anti-Gprc5B Monoclonal Antibody For Wb, Elisa - Q9Nzh0

Item Number: CM0002915

Introduction

Anti-GPRC5B Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human GPRC5B with cross-reactivity to mouse and rat. Ideal for GPCR and cell volume regulation studies.

[Learn More](#)

Parameter	Value
Product Name	GPRC5B Rabbit mAb
Remarks/Alias	RAIG2; RAIG-2; GPRC5B
Species	Human
GeneID (Human)	51704
GeneID	51704
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 304-403 of human GPRC5B (NP_057319.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NZH0
Protein Weight	45 kDa
Shipping	Ice bag