

Anti-Granulin Rabbit Monoclonal Antibody For Wb, Elisa - P28799

Item Number: CM0006381



Introduction

Rabbit monoclonal anti-Granulin (PGRN/GP88) antibody validated for WB and ELISA, cross-reactive with human and mouse. Targets lysosomal function and inflammation pathways; ideal for cancer and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Granulin Rabbit mAb
Remarks/Alias	GEP; GP88; PEPI; PGRN; CLN11; PCDGF; Granulin
Species	Human
GeneID (Human)	2896
GeneID	2896
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 494-593 of human Granulin (NP_002078.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28799
Protein Weight	64kDa
Shipping	Ice bag