

Anti-Granulin Rabbit Polyclonal Antibody For Wb And Elisa - P28799

Item Number: CM0028997



Introduction

High-quality rabbit polyclonal antibody targeting human Granulin (PGRN/GP88), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studies on lysosomal function, inflammation, and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	Granulin Rabbit pAb
Remarks/Alias	GEP; GP88; PEPI; PGRN; CLN11; PCDGF; Granulin
Species	Human
GeneID (Human)	2896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-593 of human Granulin (NP_002078.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28799
Protein Weight	64kDa
Shipping	Ice bag