

Anti-Grb14 Rabbit Polyclonal Antibody For Wb, Elisa - Q14449

Item Number: CM0025408



Introduction

GRB14 rabbit polyclonal antibody validated for Western blot and ELISA applications. Reacts with human, mouse, and rat species. Suitable for insulin signaling research and receptor interaction studies. Predicted molecular weight: 61kDa.

[Learn More](#)

Parameter	Value
Product Name	GRB14 Rabbit pAb
Remarks/Alias	GRB14
Species	Human
GeneID (Human)	2888
GeneID	2888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human GRB14 (NP_004481.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14449
Protein Weight	61kDa
Shipping	Ice bag