

## Anti-Grik2 Polyclonal Antibody For Wb, Elisa - Q13002

Item Number: CM0023503



### Introduction

Rabbit polyclonal antibody against human GRIK2 (Glutamate receptor ionotropic, kainate 2), suitable for WB and ELISA with cross-reactivity to mouse. Validated for neuroscience research on glutamate signaling and synaptic transmission.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | GRIK2 Rabbit pAb                                                                                                  |
| Remarks/Alias    | EAA4; GLR6; MRT6; GLUK6; GLUR6; GluK2; NEDLAS; GRIK2                                                              |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 2898                                                                                                              |
| GeneID           | 2898                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 30-300 of human GRIK2 (NP_068775.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q13002                                                                                                            |
| Protein Weight   | 103kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |