


Anti-Grm8 Polyclonal Antibody For Wb And Elisa - 000222

Item Number: CM0028162



Introduction

Rabbit polyclonal antibody targeting metabotropic glutamate receptor 8 (GRM8). Validated for Western blot and ELISA, reacts with human species, detects ~102 kDa protein. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	GRM8 Rabbit pAb
Remarks/Alias	GLUR8; mGlu8; GPRC1H; MGLUR8; GRM8
Species	Human
GeneID (Human)	2918
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GRM8 (NP_000836.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00222
Protein Weight	102kDa
Shipping	Ice bag