


Anti-Gsdmd (Full Length+C Terminal) Rabbit Monoclonal Antibody For Wb, Elisa - P57764

Item Number: CM0002808



Introduction

Rabbit monoclonal antibody targeting human GSDMD (full-length and C-terminal fragment). Suitable for Western blot and ELISA applications, enabling pyroptosis, inflammasome, and cell death research.

[Learn More](#)

Parameter	Value
Product Name	GSDMD (Full length+C terminal) Rabbit mAb
Remarks/Alias	DF5L, DFNA5L, FKSG10, GSDMDC1, GSDMD
Species	Human
GeneID (Human)	79792
GeneID	79792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 346-484 of human GSDMD (NP_079012.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P57764
Protein Weight	53kDa
Shipping	Ice bag