

Anti-Gsk3A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49840

Item Number: CM0016667



Introduction

Rabbit polyclonal antibody against GSK3 α , validated in WB, IF/ICC, ELISA. Reacts with human, mouse. Ideal for studying Wnt signaling, glycogen metabolism, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	GSK3 α Rabbit pAb
Remarks/Alias	GSK3 alpha; GSK3A; GSK3 α
Species	Human
Gene ID (Human)	2931
Gene ID	2931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 409-483 of human GSK3 α (NP_063937.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49840
Protein Weight	51kDa
Shipping	Ice bag