


## Anti-Gsk3B Monoclonal Antibody For Wb, Ip, If-P, Elisa - P49841

Item Number: CM0005734



### Introduction

Recombinant monoclonal GSK3 $\beta$  antibody validated for WB, IP, IF-P, and ELISA. Recognizes human, mouse, rat GSK3 $\beta$ . Suitable for studies on Wnt signaling, insulin resistance, and neurodegeneration.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] GSK3 $\beta$ Rabbit mAb                                                                                                |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 2932                                                                                                                                  |
| Gene ID          | 2932                                                                                                                                  |
| Immunogen        | Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 321-420 of human GSK3 $\beta$ (NP_001139628.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, IP, IF-P, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | P49841                                                                                                                                |
| Protein Weight   | 47kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |