


Anti-Gsta3 Rabbit Polyclonal Antibody For Wb, If, Elisa - Q16772

Item Number: CM0020456



Introduction

High-quality anti-GSTA3 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat GSTA3 (25 kDa). Involved in glutathione conjugation and steroid hormone biosynthesis. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	GSTA3 Rabbit pAb
Remarks/Alias	GTA3, GSTA3-3, GSTA3
Species	Human
GeneID (human)	2940
GeneID	2940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human GSTA3 (NP_000838.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16772
Protein Weight	25kDa
Shipping	Ice bag