

Anti-Gstk1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2Q3

Item Number: CM0020974



Introduction

CamelBio offers GSTK1 Rabbit Polyclonal Antibody, validated by knockout, for WB, IF/ICC, ELISA. Detects human GSTK1, a peroxisomal glutathione S-transferase with ~25kDa molecular weight. Ideal for oxidative stress and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GSTK1 Rabbit pAb
Remarks/Alias	GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13; GSTK1
Species	Human
GeneID (Human)	373156
GeneID	373156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human GSTK1 (NP_057001.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y2Q3
Protein Weight	25kDa
Shipping	Ice bag