

Anti-Gstm4 Polyclonal Antibody For Wb, Elisa - Q03013

Item Number: CM0026364



Introduction

Rabbit polyclonal antibody targeting human GSTM4 (glutathione S-transferase Mu 4) for WB and ELISA. Cross-reacts with mouse. Suitable for detoxification, leukotriene biosynthesis, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	GSTM4 Rabbit pAb
Remarks/Alias	GTM4; GSTM4-4; GSTM4
Species	Human
Gene ID (Human)	2948
Gene ID	2948
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human GSTM4 (NP_000841.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q03013
Protein Weight	26kDa
Shipping	Ice bag