


Anti-Gstt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30711

Item Number: CM0018586



Introduction

Rabbit polyclonal antibody against human GSTT1 for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects a 27kDa band. Ideal for oxidative stress research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	GSTT1 Rabbit pAb
Remarks/Alias	GSTT1
Species	Human
GeneID (Human)	2952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of GSTT1 (NP_000844.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30711
Protein Weight	27kDa
Shipping	Ice bag