

Anti-Gstz1 Rabbit Polyclonal Antibody For Wb, Elisa - 043708

Item Number: CM0029834



Introduction

GSTZ1 Rabbit Polyclonal Antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Detects maleylacetoacetate isomerase, a bifunctional enzyme in phenylalanine and tyrosine catabolism.

[Learn More](#)

Parameter	Value
Product Name	GSTZ1 Rabbit pAb
Remarks/Alias	MAI, MAAI, MAAID, GSTZ1-1, GSTZ1
Species	Human
Gene ID (Human)	2954
Gene ID	2954
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-216 of human GSTZ1 (NP_665877.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43708
Protein Weight	24kDa
Shipping	Ice bag