


Anti-Gtpbp4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bze4

Item Number: CM0027715



Introduction

Affinity-purified rabbit polyclonal antibody targeting GTPBP4, suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Supports ribosomal biogenesis and TP53 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GTPBP4 Rabbit pAb
Remarks/Alias	NGB; CRFG; NOG1; GTPBP4
Species	Human
Gene ID (Human)	23560
Gene ID	23560
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 385-634 of human GTPBP4 (NP_036473.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZE4
Protein Weight	74kDa
Shipping	Ice bag