

Anti-Gtpbp8 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N3Z3

Item Number: CM0018879

Introduction

Rabbit polyclonal antibody against human GTPBP8, suitable for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, and rat samples. Ideal for studying GTP-binding protein 8 function.

[Learn More](#)

Parameter	Value
Product Name	GTPBP8 Rabbit pAb
Remarks/Alias	HSPC135; GTPBP8
Species	Human
GeneID (Human)	29083
GeneID	29083
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human GTPBP8 (NP_054889.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3Z3
Protein Weight	32kDa
Shipping	Ice bag