


Anti-Gyg1 Rabbit Polyclonal Antibody For Wb And Elisa - P46976

Item Number: CM0018118



Introduction

Rabbit polyclonal antibody against human glycogenin-1 (GYG1), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for glycogen biosynthesis and glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	GYG1 Rabbit pAb
Remarks/Alias	GYG; GSD15; GYG1
Species	Human
GeneID (human)	2992
GeneID	2992
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-333 of human GYG1 (NP_001171649.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46976
Protein Weight	39kDa
Shipping	Ice bag