

Anti-Habp4 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Jvs0

Item Number: CM0020693

Introduction

Validated for WB, ELISA with cross-reactivity to Mouse, Rat. Targets intracellular hyaluronan-binding protein 4, a ribosome hibernation and transcription regulator. Ideal for cancer and transcriptional research.

[Learn More](#)

Parameter	Value
Product Name	HABP4 Rabbit pAb
Remarks/Alias	IHABP4, IHABP-4, Ki-1/57, SERBP1L, HABP4
Species	Human
Gene ID	22927
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human HABP4 (NP_055097.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5JVS0
Protein Weight	46kDa
Shipping	Ice bag