

Anti-Hapln1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P10915

Item Number: CM0030152

Introduction

Rabbit polyclonal antibody against human HAPLN1, validated in WB, IHC-P, ELISA, reactive with human, mouse, rat. Suitable for cartilage matrix research and proteoglycan aggregate studies.

[Learn More](#)

Parameter	Value
Product Name	HAPLN1 Rabbit pAb
Remarks/Alias	CRT1; CRTL1; HAPLN1
Species	Human
GeneID (human)	1404
GeneID	1404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-354 of human HAPLN1 (NP_001875.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10915
Protein Weight	40kDa
Shipping	Ice bag