

Anti-Has2 Rabbit Polyclonal Antibody For Wb, Elisa - Q92819

Item Number: CM0026165

Introduction

High-quality Rabbit Polyclonal Antibody against human Hyaluronan synthase 2 (HAS2), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for research on hyaluronan synthesis and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	HAS2 Rabbit pAb
Remarks/Alias	HAS2
Species	Human
GeneID (Human)	3037
GeneID	3037
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 67-186 of human HAS2 (NP_005319.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92819
Protein Weight	64kDa
Shipping	Ice bag