

Anti-Hcrtr1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O43613

Item Number: CM0017230



Introduction

High-affinity rabbit polyclonal antibody targeting human HCRTR1 (Orexin receptor type 1). Validated for WB, IF-P, ELISA, with cross-reactivity to mouse and rat. Suitable for neuroscience and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	HCRTR1 Rabbit pAb
Remarks/Alias	OX1R; OXR1; ORXR1; HCRTR1
Species	Human
Gene ID (Human)	3061
Gene ID	3061
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human HCRTR1 (NP_001516.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43613
Protein Weight	48kDa
Shipping	Ice bag