

Anti-Hd3A Rabbit Polyclonal Antibody For Elisa And Wb - Q93Wi9

Item Number: CM0020915



Introduction

Rabbit polyclonal antibody against rice florigen Hd3a (OsJ_20191). Validated for WB and ELISA, detects *Oryza sativa* Hd3a. Ideal for flowering time and plant developmental studies.

[Learn More](#)

Parameter	Value
Product Name	Hd3a Rabbit pAb
Remarks/Alias	Hd3a; OsJ_20191
Species	<i>Oryza sativa</i>
Gene ID	4340185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of <i>Oryza sativa</i> Hd3a. (Q93Wi9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q93Wi9
Protein Weight	20kDa
Shipping	Ice bag