

Anti-Hdac11 Rabbit Monoclonal Antibody For Wb, Elisa - Q96Db2

Item Number: CM0004123

Introduction

Rabbit monoclonal antibody targeting human HDAC11, validated for WB and ELISA. Reacts with human, mouse, and rat. Essential tool for epigenetic research, chromatin remodeling, and transcriptional regulation studies.

[Learn More](#)

Parameter	Specification
Product Name	HDAC11 Rabbit mAb
Remarks/Alias	HD11
Species	Human
Gene ID (Human)	79885
Gene ID	79885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-347 of human HDAC11 (NP_079103.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DB2
Protein Weight	39 kDa
Shipping	Ice bag