


Anti-Heatr1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H583

Item Number: CM0014389



Introduction

High-quality rabbit polyclonal antibody against human HEATR1 (BAP28/UTP10). Validated for Western blot and ELISA, this antibody is essential for exploring ribosome biogenesis, nucleolar function, and Myc-related pathways.

[Learn More](#)

Parameter	Value
Product Name	HEATR1 Rabbit pAb
Remarks/Alias	BAP28; UTP10; HEATR1
Species	Human
Gene ID	55127 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HEATR1 (NP_060542.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H583
Protein Weight	242kDa
Shipping	Ice bag