

Anti-Heparanase 1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y251

Item Number: CM0029790



Introduction

High-quality rabbit polyclonal antibody against human Heparanase 1 (HPA1), validated for WB, IF/ICC, and ELISA. Detects precursor and cleaved forms; cross-reacts with mouse and rat. Ideal for ECM degradation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Heparanase 1 Rabbit pAb
Remarks/Alias	HPA; HPA1; HPR1; HSE1; HPSE1; Heparanase 1
Species	Human
GeneID (Human)	10855
GeneID	10855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-275 of human Heparanase 1 (NP_001092010.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y251
Protein Weight	61kDa
Shipping	Ice bag