

Anti-Hk1 Polyclonal Antibody For Wb, Elisa - P19367

Item Number: CM0016757



Introduction

Anti-HK1 polyclonal antibody (pAb) for human, mouse, rat. Validated for WB and ELISA. Reacts with hexokinase-1, a central glycolysis enzyme and innate immunity receptor. Ideal for metabolism and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	HK1 Rabbit pAb
Remarks/Alias	HK; HKD; HKI; HXK1; NMSR; RP79; HMSNR; HK1-ta; HK1-tb; HK1-tc; NEDVIBA; hexokinase; HK1
Species	Human
GeneID (human)	3098
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human HK1 (NP_000179.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19367
Protein Weight	102kDa
Shipping	Ice bag