


## Anti-Hltf Rabbit Polyclonal Antibody For Wb, Elisa - Q14527

Item Number: CM0026200



### Introduction

Rabbit polyclonal antibody against human HLTF (DNA-dependent ATPase/E3 ubiquitin-protein ligase), validated for WB and ELISA. Targets a key chromatin remodeler involved in DNA damage tolerance, transcription regulation, and DNA methylation maintenance. Ideal for DNA repair and epigenetics research.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | HLTF Rabbit pAb                                                                                                                          |
| Remarks/Alias    | ZBU1; HLTf1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3; HLTf                                                                               |
| Species          | Human                                                                                                                                    |
| GeneID (Human)   | 6596                                                                                                                                     |
| GeneID           | 6596                                                                                                                                     |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 726-1009 of human HLTF (NP_620636.1) |
| Source           | Rabbit                                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                                    |
| Application      | WB, ELISA                                                                                                                                |
| Cross Reactivity | Human                                                                                                                                    |
| SWISS            | Q14527                                                                                                                                   |
| Protein Weight   | 114kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                  |