


Anti-Hmga1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17096

Item Number: CM0023704



Introduction

Rabbit polyclonal antibody against human HMGA1, validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Targets high mobility group protein HMGI/HMGI-Y involved in transcription regulation and chromatin structure.

[Learn More](#)

Parameter	Value
Product Name	HMGA1 Rabbit pAb
Remarks/Alias	HMG-R; HMGIY; HMGA1A; HMGA1
Species	Human
Gene ID (Human)	3159
Gene ID	3159
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HMGA1 (NP_665908.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17096
Protein Weight	12 kDa
Shipping	Ice bag