


Anti-Hnrnp E2/Pcbp2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q15366

Item Number: CM0003879



Introduction

Rabbit monoclonal antibody targeting human hnRNP E2/PCBP2. Validated for WB, IP, ELISA with cross-reactivity to mouse and rat. Ideal for antiviral signaling, RNA metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP E2/PCBP2 Rabbit mAb
Remarks/Alias	HNRPE2; HNRNPE2; hnRNP-E2; hnRNP E2/PCBP2
Species	Human
GeneID (Human)	5094
GeneID	5094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-350aa of human hnRNP E2/PCBP2 (NP_001122383.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15366
Protein Weight	39 kDa
Shipping	Ice bag