

Anti-Hnrrnp F Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P52597

Item Number: CM0007661



Introduction

Rabbit monoclonal antibody against human hnRNP F. Validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for studying pre-mRNA processing and alternative splicing.

[Learn More](#)

Parameter	Value
Product Name	hnRNP F Rabbit mAb
Remarks/Alias	HNRNPF; HNRPF; OK/SW-cl.23; mcs94-1; heterogeneous nuclear ribonucleoprotein F; hnRNP F
Species	Human
Gene ID	3185
Immunogen	Synthetic peptide derived from human hnRNP F
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52597
Protein Weight	46kDa
Shipping	Ice bag