

Anti-Hnrnp U Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q00839

Item Number: CM0009607

Introduction

Rabbit monoclonal antibody against human hnRNP U (Q00839) validated for WB, IP, IF/ICC, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for chromatin, RNA biology, and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP U Rabbit mAb
Remarks/Alias	SAFA; DEE54; HNRPU; SAF-A; U21.1; pp120; EIEE54; GRIP120; hnRNP U; HNRNPU-AS1
Species	Human
GeneID (Human)	3192
GeneID	3192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human hnRNP U (Q00839)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00839
Protein Weight	91 kDa
Shipping	Ice bag