

Anti-Hnrrnpl Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14866

Item Number: CM0015944



Introduction

Rabbit polyclonal antibody against human HNRNPL (hnRNP-L) for Western blot, IHC, IF/ICC, and ELISA. Validated in human, mouse, rat. Suitable for RNA splicing and neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	HNRNPL Rabbit pAb
Remarks/Alias	HNRPL; hnRNP-L; P/OKcl.14; HNRNPL
Species	Human
GeneID (Human)	3191
GeneID	3191
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14866
Protein Weight	64kDa
Shipping	Ice bag