

Anti-Homer1 Rabbit Polyclonal Antibody (Wb, If/Icc, Elisa) - Q86Ym7

Item Number: CM0023360

Introduction

Rabbit polyclonal antibody targeting human HOMER1 (SYN47), a postsynaptic scaffolding protein. Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse. Ideal for synaptic plasticity and calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HOMER1 Rabbit pAb
Remarks/Alias	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
Species	Human
GeneID (Human)	9456
GeneID	9456
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 71-170 of human HOMER1 (NP_004263.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86YM7
Protein Weight	40kDa
Shipping	Ice bag