

Anti-Hopx Polyclonal Antibody For Wb, If-P, Elisa - Q9Bpy8

Item Number: CM0024168

Introduction

Rabbit polyclonal antibody targeting HOPX (homeodomain-only protein) validated for WB, IF-P, ELISA. Detects human, mouse, rat HOPX, a regulator of cardiac development, tumor suppression, and chaperone-mediated protein folding.

[Learn More](#)

Parameter	Value
Product Name	HOPX Rabbit pAb
Remarks/Alias	HOD; HOP; OB1; LAGY; TOTO; CAMEO; NECC1; SMAP31; HOPX
Species	Human
Gene ID	84525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human HOPX (NP_631957.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPY8
Protein Weight	8kDa
Shipping	Ice bag