

Anti-Hoxa10 Polyclonal Antibody For Wb, Elisa - P31260

Item Number: CM0026645

Introduction

Rabbit polyclonal antibody targeting human HOXA10, a homeobox transcription factor involved in developmental patterning. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying gene regulation and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	HOXA10 Rabbit pAb
Remarks/Alias	PL; HOX1; HOX1H; HOX1.8; HOXA10
Species	Human
GeneID (human)	3206
GeneID	3206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 236-336 of human HOXA10 (NP_061824.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31260
Protein Weight	42kDa
Shipping	Ice bag